



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human DEK Monoclonal Antibody

单克隆抗体

产品货号: MA01950

第三版

免疫原:Fusion protein of DEK

亚型:IgG2a

克隆:88

**描述:** Dek oncogene(DEK) is first recognized in a fusion with CAN nucleoporin protein, and it's also an autoantigen in patients with polyarticular onset juvenile rheumatoid arthritis. As a site-specific binding protein, DEK may involve in transcriptional regulation and signal transduction. Analysis by agarose gel electrophoresis, DEK showed to reduce the efficiency with which chromatin is replicated. It's one of the first components of the splicing complex, which remains associated with spliced exons dependent on prior splicing of pre-mRNA.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF