

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human DHRS4 Polyclonal Antibody

多克隆抗体

产品货号: PA03526

第三版

描述: DHRS4 (Dehydrogenase/Reductase 4) is a Protein Coding gene. Among its related pathways are Signaling by Retinoic Acid and the visual cycle I (vertebrates). GO annotations related to this gene include receptor binding and oxidoreductase activity, acting on NAD(P)H, quinone or similar compound as acceptor. An important paralog of this gene is DHRS4L2.

抗原: Synthesized peptide derived from the C-terminal region of human DHRS4

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Isoform 1 is predominantly expressed in normal cervix (at protein level). Isoform 4 is expressed in some neoplastic cervical tissues, but not in normal cervix (at protein level). Isoform 5 and isoform 6 are expressed in a few neoplastic cervical tissues.