

Aminoacyl tRNA synthetase complex-interacting multifunctional protein 2 (AIMP2)

Monoclonal Antibody

Cat.NO.: MA01075

3th Edition

Immunogen:Denatured full length native recombinant human AIMP2.

Isotype:IgG

Clone:

Description:

Form:Anti-AIMP2 antibody containing 0.05% sodium azide for preservative.

Concentration:Please see datasheet.

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF