

BPGM Anti-Human Antibody (Clone AT2E11)

Monoclonal Antibody

Cat.NO.: MA01121

3th Edition

Immunogen: Recombinant human BPGM (1-259aa) purified from E. coli.

Isotype: Mouse IgG2a heavy chain and kappa light chain

Clone: Anti-human BPGM mAb, clone AT2E11, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human BPGM protein.

Description:

Form: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol

Concentration: 1 mg/ml

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF