

**c-MPL/TPOR (AMM2) Anti-Mouse Rat IgG MoAb Biotin**  
**Monoclonal Antibody**  
**Cat.NO.: MA01130**

---

3th Edition

**Immunogen:** Recombinant protein of mouse c-MPL extracellular domain.

**Isotype:** IgG1

**Clone:** AMM2

**Description:**

**Form:** Lyophilized product in PBS containing 1 % BSA and 0.05 % NaN<sub>3</sub>

**Concentration:** Please see datasheet.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF