

Anti-CD252/OX40L (8F4) Mouse IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01155

3th Edition

Immunogen: Lysate of HTLV-1 persistently-infected cell line (MT-2)

Isotype: IgG1

Clone: 8F4

Description:

Form: Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN₃

Concentration: Please see datasheet.

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF