

**Cyclin H Antibody (Clone AT3G6)**

**Monoclonal Antibody**

**Cat.NO.: MA01214**

---

3th Edition

**Immunogen:** Recombinant human Cyclin H (1-323aa) purified from E.coli

**Isotype:** Mouse IgG2b, k

**Clone:** Anti-human Cyclin H mAb, clone AT3G6, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human Cyclin H protein.

**Description:**

**Form:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**Concentration:** 1 mg/ml

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF