

**EIF2S1 Anti-Human Antibody (Clone AT5E10)**

**Monoclonal Antibody**

**Cat.NO.: MA01232**

---

3th Edition

**Immunogen:** Recombinant human EIF2S1 (1-315aa) purified from E. coli

**Isotype:** IgG2b heavy chain and kappa light chain.

**Clone:** 5E10.

**Description:**

**Form:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**Concentration:** 1 mg/ml

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF