

**EPCAM Antibody (Clone AT37F9)**

**Monoclonal Antibody**

**Cat.NO.: MA01237**

---

3th Edition

**Immunogen:** Recombinant human EPCAM (24-265aa) purified from E. coli

**Isotype:** IgG1,k

**Clone:** Anti-human EPCAM mAb, clone AT37F9, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human EPCAM protein.

**Description:**

**Form:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**Concentration:** 1mg/ml

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF