

Gd-IgA1 (KM55) Anti-Human Rat IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01276

3th Edition

Immunogen: Human IgA1 hinge region peptide with GalNAc H-C223PST*PPT*PS*PS*TPPT*PSPS240-NH2 (*with GalNAc)

Isotype: Rat IgG2b

Clone: KM55

Description:

Form: Lyophilized product from PBS containing 1% BSA and 0.05 % NaN₃

Concentration: 1 mL deionized water will be added to the product, then its concentration comes to 100 µg/mL.

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF