

GRO/CINC-1 (6G2) (C) Anti-Rat Mouse IgM MoAb

Monoclonal Antibody

Cat.NO.: MA01293

3th Edition

Immunogen: Synthetic peptide to C-terminal of Rat GRO/CINC-1

Isotype: IgM

Clone: 6G2

Description:

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

Concentration: 500 ug/mL

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF