

Tau (Phosphorylated) (C5) Anti-Human Mouse IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01318

3th Edition

Immunogen:PHF (paired helical filament) obtained from AD brain.

Isotype:IgG1

Clone:C5

Description:

Form:Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃

Concentration:

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF