

Tenascin-C (Domain B) (4C8MS) Anti-Human Mouse IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01327

3th Edition

Immunogen: Human Recombinant Tenascin-C

Isotype: IgG1

Clone: 4C8MS

Description:

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

Concentration: Please see datasheet.

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF