

Anti- α , β -Synuclein Antibody(Clone 3B6)

Monoclonal Antibody

Cat.NO.: MA01381

3th Edition

Immunogen:Recombinant human alpha-synuclein (119-140aa) purified from E. coli

Isotype:Mouse IgG1 heavy chain and kappa light chain

Clone:3B6

Description:

Form:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

Concentration:1mg/ml

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF