

Hexokinase1 Anti-Human Antibody (Clone 4D7)

Monoclonal Antibody

Cat.NO.: MA01411

3th Edition

Immunogen: Recombinant human Hexokinase1 (1-917aa) purified from E. coli

Isotype: Mouse IgG2a heavy chain and kappa light chain

Clone: Anti-human Hexokinase mAb, clone 4D7, is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from BALB/c mice

Description:

Form: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

Concentration: 1 mg/ml

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF