

**IDH1 R132S (SMab-1) Anti-Human Mouse IgG MoAb**

**Monoclonal Antibody**

**Cat.NO.: MA01433**

---

3th Edition

**Immunogen:** Synthetic peptide of a part of human IDH1 R132S

**Isotype:** IgG1

**Clone:** SMab-1

**Description:**

**Form:** Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN<sub>3</sub>

**Concentration:** Please see datasheet.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF