



**MBD3 (Methyl-CpG-Binding Protein 3) (3A3) Anti-Human Mouse IgG MoAb**

**Monoclonal Antibody**

**Cat.NO.: MA01509**

---

3th Edition

**Immunogen:** Synthetic peptide of part of C terminal of human MBD3 (EELARDGEAPLDKA)

**Isotype:** IgG1

**Clone:** 3A3

**Description:**

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**Concentration:** Please see datasheet.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF