

Reticulocalbin-1 (TMU-6A1) Anti-Human Mouse IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01545

3th Edition

Immunogen: Recombinant protein of the C-terminal site (90 a.a.) of human Reticulocalbin-1

Isotype: IgG1

Clone: TMU-6A1

Description:

Form: Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃

Concentration:

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF