

**Anti-Neurofilament, 160K, 210K, (N21-8-4) Mouse Ascites MoAb**

**Monoclonal Antibody**

**Cat.NO.: MA01564**

---

3th Edition

**Immunogen:**Bovine Brain Neurofilament

**Isotype:**IgG1

**Clone:**N-21-8-4

**Description:**

**Form:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**Concentration:**Please see datasheet.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:**The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:**ELISA, WB, ICC/IF