

**Anti-Smad3L (S213 Phosphorylated) (5A11) Mouse IgG MoAb**

**Monoclonal Antibody**

**Cat.NO.: MA01592**

---

3th Edition

**Immunogen:** Synthetic peptide of phosphorylated Smad3L (S213).

**Isotype:** IgG1

**Clone:** 5A11

**Description:**

**Form:** Lyophilized product in PBS containing 1 % BSA and 0.05 % NaN3

**Concentration:** Please see datasheet.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF