

**SP-A (PE10) Anti-Human Mouse IgG MoAb**

**Monoclonal Antibody**

**Cat.NO.: MA01606**

---

3th Edition

**Immunogen:**Purified surfactant apoprotein A

**Isotype:**IgG2b

**Clone:**PE10

**Description:**

**Form:**Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN<sub>3</sub>

**Concentration:**Please see datasheet.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:**The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:**ELISA, WB, ICC/IF