



ORP150 (Oxygen Regulated Protein, p150) (2F07) Anti-Human Mouse IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01613

3th Edition

Immunogen:Recombinant Protein (GST-ORP150 (Leu 508 - Leu 999))

Isotype:IgG2a

Clone:2F07

Description:

Form:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

Concentration:Please see datasheet.

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF