

Parkin (1A1) Anti-Human Mouse IgG MoAb
Monoclonal Antibody
Cat.NO.: MA01642

3th Edition

Immunogen: Synthetic peptide of part of C terminal of human parkin (AYRVDERAAEQARWEAA)

Isotype: IgG1

Clone: 1A1

Description:

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

Concentration: Please see datasheet.

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF