

**Pcna Anti-Human Antibody (Clone K1A3)**

**Monoclonal Antibody**

**Cat.NO.: MA01655**

---

3th Edition

**Immunogen:**Recombinant human PCNA (1-261aa) purified from E. coli

**Isotype:**Mouse IgG3 heavy chain and kappa light chain

**Clone:**k1A3

**Description:**

**Form:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**Concentration:**1mg/ml

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:**The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:**ELISA, WB, ICC/IF