

Anti-Poly (ADP-ribose) (10H) Mouse IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01707

3th Edition

Immunogen:Purified Poly (ADP-ribose)

Isotype:IgG3, ?

Clone:10H

Description:

Form:Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN₃

Concentration:Please see datasheet.

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF