



Anti-Human/Mouse/Rat Mono-Methyl-Histone H3 (Lys79) Monoclonal Antibody

Monoclonal Antibody

Cat.NO.: MA01809

3th Edition

Immunogen:Synthetic Peptide

Isotype:IgG

Clone:3L4

Description:

Form:

Concentration:1 mg/mL

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF