



**EYFP Monoclonal Antibody**  
**Monoclonal Antibody**  
**Cat.NO.: MA01832**

---

3th Edition

**Immunogen:** Recombinant Protein

**Isotype:** IgG

**Clone:** 3B7

**Description:**

**Form:**

**Concentration:** 1mg/mL

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF