

**Anti-Human/Mouse/Rat CASP9 Monoclonal Antibody**

**Monoclonal Antibody**

**Cat.NO.: MA01886**

---

3th Edition

**Immunogen:** Fusion protein of Caspase 9/p35/p10

**Isotype:** IgG2b

**Clone:** 34

**Description:**

**Form:**

**Concentration:** None

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF