

Anti-Human CD146 Monoclonal Antibody
Monoclonal Antibody
Cat.NO.: MA01894

3th Edition

Immunogen: Fusion protein of CD146/MCAM

Isotype: IgG1

Clone: 68

Description:

Form:

Concentration: None

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF