

Anti-Human CFB Monoclonal Antibody
Monoclonal Antibody
Cat.NO.: MA01928

3th Edition

Immunogen: Fusion protein of Complement factor B

Isotype: IgM

Clone: 12

Description:

Form:

Concentration: None

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF