

**Anti-Human/Mouse/Rat/Monkey/Cow/Dog/Hamster/Pig/Rabbit/Sheep/Insect/Yeast GAPDH Monoclonal Antibody**

**Monoclonal Antibody**

**Cat.NO.: MA02343**

---

3th Edition

**Immunogen:**Synthetic Peptide

**Isotype:**IgG

**Clone:**9S5

**Description:**

**Form:**

**Concentration:**1mg/mL

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:**The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:**ELISA, WB, ICC/IF