



TAP Tag Monoclonal Antibody
Monoclonal Antibody
Cat.NO.: MA02359

3th Edition

Immunogen: Recombinant Protein

Isotype: IgG

Clone: 6B3

Description:

Form:

Concentration: 1mg/mL

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF