

Anti-Human/Mouse/Rat/Monkey/Chicken/Dog/Hamster/Rabbit beta actin Monoclonal Antibody

Monoclonal Antibody

Cat.NO.: MA02363

3th Edition

Immunogen:Synthetic Peptide

Isotype:IgG

Clone:2K4

Description:

Form:

Concentration:1mg/mL

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF