

**Anti-Mouse LE/AF Purified Anti-Mouse H-2 Monoclonal Antibody**

**Monoclonal Antibody**

**Cat.NO.: MA02418**

---

3th Edition

**Immunogen:**None

**Isotype:**Rat IgG2a

**Clone:**M1/42.3.9.8

**Description:**

**Form:**

**Concentration:**None

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:**The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:**ELISA, WB, ICC/IF