

Anti-Human/Mouse/Rat ARHGDI Monoclonal Antibody

Monoclonal Antibody

Cat.NO.: MA02434

3th Edition

Immunogen: Fusion protein of ARHGDI/aGDI

Isotype: IgG2a

Clone: 1B1G9

Description:

Form:

Concentration: None

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF