

**Anti-Human/Mouse/Rat CALB1 Monoclonal Antibody**

**Monoclonal Antibody**

**Cat.NO.: MA02441**

---

3th Edition

**Immunogen:**Fusion Protein of Calbindin-D28k

**Isotype:**IgG1

**Clone:**1F8B9

**Description:**

**Form:**

**Concentration:**None

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:**The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:**ELISA, WB, ICC/IF