

**Anti-Human/Mouse/Rat MAP2K2 Monoclonal Antibody**

**Monoclonal Antibody**

**Cat.NO.: MA02655**

---

3th Edition

**Immunogen:** Synthesized peptide derived from human MEK-2 around the non-phosphorylation site of T394.

**Isotype:** IgG

**Clone:**

**Description:**

**Form:**

**Concentration:** 1mg/mL

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF