



**Anti-Human/Mouse/Rat Acetyl-Histone H3 (Lys9) Monoclonal Antibody**

**Monoclonal Antibody**

**Cat.NO.: MA02708**

---

3th Edition

**Immunogen:** Synthetic Peptide of Acetyl-Histone H3 (K9)

**Isotype:** IgG

**Clone:** 2E7

**Description:**

**Form:**

**Concentration:** 1mg/mL

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF