

**Adiponectin Receptor 1 Anti-Human Rabbit IgG Affinity Purify
Polyclonal Antibody**

Cat.NO.: PA01010

3th Edition

Description:Adiponectin has been shown to be one of the adipokines possessing anti-diabetic, anti-atherosclerotic and anti-inflammatory actions. It was recently shown that the blood adiponectin levels were reduced in cases of obesity.

Antigen:Synthetic peptide of a part of human Adiponectin Receptor 1.

Form:Lyophilized product from 1 % BSA in PBS containing 0.05% NaN₃.

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Not cross-react with Adiponectin Receptor 2. Cross-reacts with Mouse.