

14-3-3 ? Protein Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01018

3th Edition

Description:The 14-3-3 proteins are a family of conserved regulatory molecules expressed in all eukaryotic cells. A striking feature of the 14-3-3 proteins is their ability to bind a multitude of functionally diverse signaling proteins, including kinases, phosphatases, and transmembrane receptors.

Antigen:Synthetic peptide for C-terminal of Human 14-3-3 ?.

Form:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Human 14-3-3 epsilon specific, Non-crossreact to Human 14-3-3 ?, ?, ?, ?, ?.