

**APP (597) Anti-Mouse/Rat Rabbit IgG Affinity Purify
Polyclonal Antibody**

Cat.NO.: PA01056

3th Edition

Description: Amyloid precursor protein (APP) is precursor protein of Amyloid ? which is major constituent of senile plaque in Alzheimer's disease. It is known that there are three major isoforms.

Antigen: Synthetic peptide in portion of C terminus of Mouse/Rat sAPP? (DAEFGHDSGFVVRHQK).

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Not cross-react with Human. Not cross-react with sAPP?.