

**Anti-BACE1 (42) Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01065**

---

3th Edition

**Description:** The current understanding of the role of amyloid in Alzheimer's disease (AD) has been established by a remarkable congruence of multiple disciplines: neuropathology, biochemistry, molecular biology, and epidemiological genetics.

**Antigen:** Synthetic peptide of a part of human BACE1 (RLPRETDEEPEEPGRR)

**Form:** Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Reacts with human, mouse and rat BACE1.