

Anti-BACE1 (42) Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01065

3th Edition

Description:The current understanding of the role of amyloid in Alzheimer's disease (AD) has been established by a remarkable congruence of multiple disciplines: neuropathology, biochemistry, molecular biology, and epidemiological genetics.

Antigen:Synthetic peptide of a part of human BACE1 (RLPRETDEEPEEPGRR)

Form:Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Reacts with human, mouse and rat BACE1.