

c-Met ? Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01078

3th Edition

Description: Overexpression of the hepatocyte growth factor receptor (c-Met/HGF receptor), a transmembrane tyrosine kinase encoded by the MET proto-oncogene, is involved in transformation and invasive behavior of human carcinomas and sarcomas.

Antigen: Synthetic peptide for C-terminal of Human c-Met ?

Form: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: React to c-Met beta chain.