

Dock7 Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01132

3th Edition

Description: Small GTPase of Rho family is one of signaling factors which control alteration of cellular morphology, cell division and differentiation, and plays an important role in regulation of generation and organogenesis at an individual level. Its functions are strictly regulated and its dysfunction causes many diseases such as immunological disease or neurodegenerative disease.

Antigen: Synthetic peptide of a part of Human Dock7 (KELFALHPSPDEEE)

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Cross-reacts with rat.