

Anti-Drebrin A (DAS2) Rabbit IgG Affinity Purity

Polyclonal Antibody

Cat.NO.: PA01136

3th Edition

Description: Drebrin (developmentally regulated brain protein), isolated and identified from chicken brain as a protein associated with brain development, is conserved over species, and there are two types of isoforms, drebrin A and drebrin E in mammalian such as rat and human.

Antigen: Synthetic peptide of the part of Drebrin A (FIKASDSGPSSS)

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Reacts with rat and mouse Drebrin A.