

Anti-GFAT1 (663) Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01165

3th Edition

Description:GFAT1 (Glutamine : fructose-6-phosphate amidotransferase 1) is one of nutrientsensitive phosphoprotein which was identified as a protein phosphorylated in glucose-deprived cells.

Antigen:Synthetic peptide of the C terminal part of GFAT1 (RGYDVDFPRNLAKSVTVE)

Form:Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃.

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Cross reacts with mouse and rat.