

**Anti-GFAT1 (663) Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01165**

---

3th Edition

**Description:**GFAT1 (Glutamine : fructose-6-phosphate amidotransferase 1) is one of nutrientsensitive phosphoprotein which was identified as a protein phosphorylated in glucose-deprived cells.

**Antigen:**Synthetic peptide of the C terminal part of GFAT1 (RGYDVDFPRNLAKSVTVE)

**Form:**Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN<sub>3</sub>.

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Cross reacts with mouse and rat.