

**GRO/CINC-2? (C) Anti-Rat Rabbit IgG Affinity Purify
Polyclonal Antibody**

Cat.NO.: PA01183

3th Edition

Description: Cytokine-induced neutrophil chemo attractant 1 (CINC-1) was originally purified from media conditioned by IL-1 ? stimulated rat kidney epithelioid cells (NRK-52E).

Antigen: Synthetic peptide for C-Terminal of rat GRO/CINC-2? conjugating bovine thyroglobulin.

Form: Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Rat GRO/CINC-2 ? specific. Non-cross react with Rat GRO/CINC-1, Rat GRO/CINC-2? and Rat GRO/CINC-3