

**GRO β /MGSA (C) Anti-Human Rabbit IgG Affinity Purify
Polyclonal Antibody**

Cat.NO.: PA01188

3th Edition

Description: GRO β /MGSA as IL-8 family has the high homology with IL-8 in amino acid sequences. It is reported that it also has the similar configuration in 3D structural analysis and the almost same level of leukocyte chemotaxis with IL-8.

Antigen: Synthetic peptide in portion of C-terminus of Human GRO β /MGSA.

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1 μ g/ml.

Specificity: Non-Crossreacts to Human GRO β , Human IL-8, Human IP-10, Rat GRO/CINC-1, Rat GRO/CINC-2 β , Rat GRO/CINC-2 α and Rat GRO/CINC-3.